



C.S.A. SERVICE Srl - Via Farini, 8
20030 SENAGO - MI - Italy
Tel. (+39) 02 99813304
Fax (+39) 02 99815209
<http://www.csamedical.it/>
e-mail: info@csamedical.it

DISPOSABLE BIB WITH POCKET

BAV01 - DISPOSABLE BIB WITH POCKET

In accordance with EU directive 93/42/EEC and its modifications

Class: 1°

Type: disposable not sterile.

Characteristics

Perfect for body cleaning of old and disabled people. Very absorbent and totally water resistant.

It can be used in both sides. Excellent resistance to longitudinal and transversal traction.

Provided with pocket.

Weight: 11,4 g. each bib.

Colour: light blue.

Material

Absorbent side in tissue: 24 g/m².

Water resistant side in blue polyethylene 0,014 mm thick.

Adhesive: water based.

Measures: 37 x 66 cm.

Confections and packing

125 pcs / PE-bag, 500 pcs / box, 54 boxes / Eur-pallet.

