



C.S.A. SERVICE Srl - Via Farini, 8
20030 SENAGO - MI - Italy
Tel. (+39) 02 99813304
Fax (+39) 02 99815209
<http://www.csamedical.it/>
e-mail: info@csamedical.it

WASHING GLOVE WITH PE

MAN02 - WASHING GLOVE WITH PE

In accordance with EU directive 93/42/EEC and its modifications



Class: 1°

Type: disposable not sterile.

Characteristics

Perfect for body cleaning of old and disabled people. Totally water resistant.

Excellent resistance to longitudinal and transversal traction.

Weight: 5,5 g. each washglove. Colour: light blue.

Material

External in Airlaid paper: 50 g/m². Internal in polyethylene: blu. Adhesive: water based.

Measures: 16 x 24 cm.

Confections and packing

100 pcs / PE-bag, 1000 pcs / box, 50 boxes / Eur-pallet.